

KARMI Blood Bags

Single and Multiple
Blood Collection System

 **KAWASUMI LABORATORIES, INC.**

Visit on our Web site at www.kawasumi.jp

KAWASUMI GROUP

MANUFACTURER

KAWASUMI LABORATORIES, INC.

Shinagawa Intercity B, 9th Floor
2-15-2, Konan, Minato-ku, Tokyo
ZIP: 108-6109 Japan
TEL: +81-3-5769-2605 Fax: +81-3-5769-276

KAWASUMI LABORATORIES AMERICA, INC.

P.O. Box 24355 Tampa, Florida U.S.A. 33623-4355

Kawasumi Laboratories (Thailand) Co., Ltd. (KLT-3B)

Nava Nakorn Industrial Promotion Zone, 55/26 MU 13, Phanon Yothin Rd., KM-46, Tambon Khlong Nueng,
Amphoe Khlong Luang, Changwat Pathum Thani 12120, Thailand

NAMSIN TRADING CO., LTD.

9th Floor Maneeya Center Building 518/5 Ploenchit Road, Lumpini, Patumwan, Bangkok 10330, Thailand

Manufacturing Plant

Saiki Plant : 1077, Koda, Yayoi, Saiki, Oita 876-0121, Japan

Navanakorn Plant : Nava Nakorn Industrial Promotion Zone, 55/26 MU 13, Phanon Yothin Rd.,
KM-46, Tambon Khlong Nueng, Amphoe Khlong Luang,
Changwat Pathum Thani 12120, Thailand

EC REGULATORY REPRESENTATIVE

MPS Medical Product Service GmbH

Borngasse 20, 35619 Braunfels, Germany

CE 0123

KARMI Blood Bag

2013-04

2010-05-20

CE 0123

KAWASUMI LABORATORIES, INC.

1005207AS

Distributor :

KARMI® BLOOD BAGS

SINGLE BAG



Type	Collection Bag Capacity	Donor Needle / Anticoagulant	Packaging
Single KBS-250 CA7 KBS-350 CA7 KBS-450 CA7	250 ml 350 ml 450 ml	17 G/ CPDA-1	- Box 100 pcs - 20 aluminium pouch @ 5 bags in individual sterile plastic pouch

DOUBLE BAG



Type	Collection Bag Capacity	Donor Needle / Anticoagulant	Packaging
Double KBD-250 CA7 KBD-350 CA7 KBD-450 CA7	250 ml 350 ml 450 ml	17 G/ CPDA-1	- Box 60 Pcs - 12 aluminium pouch @ 5 bags in individual sterile plastic pouch

TRIPLE BAG



Type	Collection Bag Capacity	Donor Needle / Anticoagulant	Packaging
Triple KBT-350 CA7 KBT-450 CA7	350 ml 450 ml	17 G/ CPDA-1	- Box 50 Pcs - 10 aluminium pouch @ 5 bags in individual sterile plastic pouch

Specifications :

- Individual packing in sterile plastic pouch
- 17 g hypodermic needle
- CPDA-1 anticoagulant
- 5 days platelets storage bag in triple bags
- Ultra sharp needle
- Snap tip port
- Extremely elastic tubing
- Detachable port protector
- Ultra smooth inner bag material to minimize cellular components damage
- Rounded corner assists smooth transfer of blood during separation and transfusion