



DiaCent-12 Bloodbank Centrifuge

Product Specifications

- Removable centrifuge head for easy cleaning and disinfection
- Time counter displayed in minutes:seconds
- Swing out rotor 45 accomodating both good sample separation and reading agglutination
- Centrifuge head for 12 tubes 10 x 75 mm or 12 x 75 mm
- Monitored by microprocessor
- Very quiet
- Rpm, time and functions are displayed (LCD) in English, French, German, Italian, Spanish and Portuguese
- Adjustable parameters or 10 different programs
- Max. RCF 1,006 g and 3,000 rpm
- Dimensions: 26 x 32 x 18 cm (W x D x H)
- Weight: 4.8 kg
- 110-240 V/50-60 Hz
- CE compliant, according to IVD Directive 98/79/EC



DiaCent-12 Centrifuge

Ordering Information

Catalog No.	Description
840000	DiaCent-12 (complete with head Catalog No. 840010)1 unit
840010	Accessory head for 12 tubes (10-12 mm diameter)1 unit

The Complete Solution for Safe Transfusion

