

ID - Centrifuge 12 S II

ID-Centrifuge 12 S II



ID-Centrifuge 12 S II

Product specifications

- Max capacity for 12 ID-cards
- Monitored by microprocessor
- Automatic balance control of the centrifuge head
- Automatic error message generation
- Removable centrifuge head
- Pre-fixed centrifugation speed and time at 1030 rpm \pm 5 rpm (85 g \pm 1 g) for 10 minutes
- Rpm or G-force and time are displayed by LCD in English, French, German, Italian, Spanish or Portuguese
- Applicable for all immunohaematological tests using DiaMed ID cards, including those using AHG reagent
- Dimension: 30 x 36 x 18 cm (W x D x H)
- Weight: 6.2 kg
- Automatic voltage 110V - 240V / 50 Hz - 60 Hz
- CE compliant, according to IVD Directive 98 / 79 / EC

Ordering Information

Catalog No.	Description
009550V	ID-Centrifuge 12 S II (complete with head)1 unit.

The Complete Solution for Safe Transfusion