



ID-Diluent 2

ID-Diluent 2 is a modified low ionic strength solution made for the ID-Micro Typing System, for preparing 5% red cell suspensions for blood grouping as well as 0.8%-1% red cell suspensions for crossmatch, autocontrol, direct anti-human globulin test, blood grouping of newborns and test cells prepared in the laboratory.



ID-Diluent 2

Product specifications

- Modified LISS (Low Ionic Strength Solution) for red blood cell suspensions
- Inhibits non-specific reactions due to complement activation
- Suitable for use with anticoagulated or native blood samples
- Storage : 2 - 8°C

Ordering Information

Catalog No.	Description
009260	Modified LISS for red cell suspensions 2 x 100 ml
009260	Modified LISS for red cell suspensions 1 x 500 ml

The Complete Solution for Safe Transfusion