



# ID-Dispenser

## Product specifications

- For serial dispensing of 0.5 ml or 1 ml Diluent 1 and Diluent 2
- Volume manually adjustable for dispensing 0.5 ml or 1 ml
- Interchangeable colour code ( red and green ) for Diluent 1 and Diluent 2



ID-Dispenser

## Ordering Information

| Catalog No. | Description        |             |
|-------------|--------------------|-------------|
| 009612V     | ID-Dispenser Red   | .....1 unit |
| 009614V     | ID-Dispenser Green | .....1 unit |

## The Complete Solution for Safe Transfusion



*Bio-Rad Laboratories is a multinational manufacturer and distributor of clinical diagnostics and life science research products.*