


ID - Incubator 37 S II



# ID-Incubator 37 S II



ID-Incubator 37 S II

## Product specifications

- For the incubation of ID-Cards at 37°C
- Max capacity of 36 ID-Cards or 1 ID-Centrifuge head ( 6, 12 or 24 ID-Cards )
- Monitored by microprocessor
- Automatic error message generation
- Standard program: 15 minutes at 37°C ( ± 2°C )
- 3 individual timers for 3 separate batches
- Programmable incubation time from 1 - 99 minutes
- Time and temperature are displayed in LCD in English, French, German, Italian, Spanish or Portuguese
- Applicable for all immunohematological tests using DiaMed ID-Cards
- Dimension: 36 x 42 x 21 cm ( W x D x H )
- Weight: 9 kg
- Automatic voltage 110V - 240V / 50 Hz - 60 Hz
- CE compliant, according to IVD Directive 98/79/EC

## Ordering Information

| Catalog No. | Description                      |
|-------------|----------------------------------|
| 009691      | ID-Incubator 37 S II .....1 unit |

The Complete Solution for Safe Transfusion



Bio-Rad Laboratories is a multinational manufacturer and distributor of clinical diagnostics and life science research products.