

Platelet Agitators

Innovation can be beautiful.

PF48i



PF96i



Benefits

Delrin acetal resin rollers and glides for smooth agitation

The added safety of Bacteria-Resistant Powder Coating

Perforated drawers with non-slip texture provide uniform air circulation for your platelets

Built-in motion alarm

Features

- Made from galvanized steel with bacteria-resistant powder coating
- Superior engineered frame for maximum rigidity
- One-piece, removable perforated drawers with non-slip texture provide uniform air circulation for your platelets
- Pull-out drawers for easy access to platelets
- Built-in motion alarm
- Single fan for forced air cooling of motor compartment
- Label holders provided for each agitator drawer accommodate magnetic and adhesive labels

SPECIFICATIONS	PF48i	PF96i
Material	Galvanized steel with bacteria-resistant powder coating	
Alarm	Built-in motion alarm	
Drawers	One-piece, removable, perforated, texture coated non-slip drawers	
	7 pull-out, 1 stationary	7 pull-out, 1 stationary
Capacity	48 random bags / 16 apheresis bags	96 random bags / 32 apheresis bags
Dimensions (wxhxd) (mm)	458 x 388 x 369	851 x 388 x 369
Platelet incubator compatibility	PC900	PC1200
Electrical	230V 50/60Hz 0.2 FLA	230V 50/60Hz 0.2 FLA
Net Weight	23 kg	37 kg

Certified to applicable UL and CSA standards by a NRTL
 Meets the standards established by the AABB, ARC and FDA for blood storage

